

*It's for everyone.*



K-0391 Black



K-0391-G White



K-0391-V Violet

Inami Universal Trial Frame K-0391 is designed to maximise examinee's comfort by minimising the frame weight. This highly durable instrument features quick and easy adjustment for both PD and corneal alignment, while providing a sure mounting for trial lenses.



## K-0391 Series

Inami Universal Light Weight Trial Frame

**PD** : 48 ~ 80mm ( $\pm 0.8$ mm)

**Cornea measuring plate** : 2mm interval

**Frame adjustable angles** : 35° (upward 30°  $\pm 5^\circ$ , downward 5° + 5° - 2°)

**Frame adjustable length** : 35mm

**Nose pad adjustable angles** : 360°

**Nose pad adjustable length** : 22mm  $\pm 1.5$ mm

**Maximum available lenses** : 10 pcs

**Axis of lens** : 0°  $\pm 5^\circ$ , 180°  $\pm 5^\circ$

**Available colour** : Black (K-0391), Violet (K-0391-V), White (K-0391-G)