


OPHTHALMIC INSTRUMENTS

INAMI GENERAL PRODUCT CATALOGUE



SLIT LAMP MICROSCOPE



L-0185 LED / Halogen

Slit lamp microscope
Greenough



L-0187 LED / Halogen

Slit lamp microscope
3 step



L-0189 LED / Halogen

Slit lamp microscope
5 step



L-0240 LED / Halogen

Slit lamp microscope
Zoom



L-0229-E LED / Halogen

Digital Camera Compatible
Slit lamp "EGG", 5 step

* Colour selectable base sheet (Blue / pink) is included in a standard set of "PHAROS" series.

	Standard / Option	L-0185	L-0187	L-0189	L-0240	L-0229-E
Type		Greenough Type Stereoscopic Microscope	Convergent Type Binocular Microscope	Convergent Type Binocular Microscope	Convergent Type Binocular Microscope	Convergent Type Binocular Microscope
Magnification Changer		2 step by changing knob	3 step by drum rotation	5 step by drum rotation	Continuous Zoom	5 step by drum rotation
Stereo angle		13°	12°	12°	12°	12°
Magnification Ratio (Field of View)		10x (20.13mm) 16x (16.13mm)	10x (28.5mm) 16x (14mm) 25x (7mm)	6x (38mm) 10x (24mm) 16x (15mm) 25x (9mm) 40x (6mm)	Continuous from 10x ~ 30x (22mm~8mm)	6x (38mm) 10x (24mm) 16x (15mm) 25x (9mm) 40x (6mm)
Slit Width		0-10 mm	0-10 mm	0-10 mm	0-10 mm	0-14 mm
Filters	Standard	Cobalt blue, Heat absorption, 1/2 ND (Halogen), Incandescent (LED), Red-free, Yellow (Observation path)	Cobalt blue, Heat absorption, 1/2 ND (Halogen), Incandescent (LED), Red-free	Cobalt blue, Heat absorption, 1/2 ND (Halogen), Incandescent (LED), Red-free	Cobalt blue, Heat absorption, 1/2 ND (Halogen), Incandescent (LED), Red-free	Cobalt blue, Heat absorption, 1/2 ND (Halogen), Incandescent (LED), Red-free
	Option	—	YF-57 adaptor built-in yellow filter	YF-57 adaptor built-in yellow filter	YF-57 adaptor built-in yellow filter	YF-57 adaptor built-in yellow filter
Colour selectable base		Selectable	Selectable	Selectable	Selectable	—
Weight		18.5 Kgs	18.5 Kgs	18.5 Kgs	18.5 Kgs	23 Kgs
Dimensions		LED: 550(W)x400(D)660(H)mm Halogen: 550(W)x400(D)690(H)mm	LED: 550(W)x400(D)660(H)mm Halogen: 550(W)x400(D)690(H)mm	LED: 550(W)x400(D)660(H)mm Halogen: 550(W)x400(D)690(H)mm	LED: 550(W)x400(D)660(H)mm Halogen: 550(W)x400(D)690(H)mm	LED: 550(W)x400(D)660(H)mm Halogen: 550(W)x400(D)690(H)mm

OPTIONAL ACCESSORY



L-5110

Applanation tonometer
swing arm type



L-0541DC-5

Digital video camera
beam splitter 5M with
filing software "Pixel Filer"

* L-0541DC-12(12M) is also available



L-0542-CB

Co-axial observation tube
(Binoculars)



L-0541

Beam splitter device

OPERATING MICROSCOPE



L-0955 / L-0940 series

Portable operation microscope
with carrying case



L-0955SD LED

Portable operation microscope
with carrying case



L-0955XD LED

Portable operation microscope
with carrying case



L-0940SD Halogen

Portable operation microscope
with carrying case



K-1422

Optional pole and base for
L-0940 / K-0955 series



L-0995AZ

LED Operation microscope with electric focusing and X-Y device



L-0990A

Operation microscope with electric focusing and X-Y device



L-0990

Operation microscope with electric focusing



L-0970

Multi operation microscope with electric focusing



Smooth vertical movement



Flexible movement of microscope head


L-0955SD	LED	Zoom(5x-15x) Objective lens:F175mm
L-0955XD	LED	Zoom(5x-15x),Motorised focus Objective lens:F175mm
L-0940SD	Halogen	Zoom(6x-18x) Objective lens:F175mm
L-0995AZ	LED	Zoom(4.4x-26.4x),Motorised focus Objective lens:F200mm
L-0990A	Halogen	5-steps(4x,6x,10x,16x,25x),X-Y device Objective lens:F175mm
L-0990	Halogen	5-steps(4x,6x,10x,16x,25x) Objective lens:F175mm
L-0970	Halogen	5-steps(4x,6x,10x,16x,25x),X-Y device Objective lens:F250mm for ENT(Multi purpose)

Optional accessories:

Beam Splitter, CCD Camera Adaptor, Assistant Microscope, Observation Tube, Oblique Illuminator, 45° Adaptor (L-0970)
Built-in C-mount for camera connection (L-0940 / K-0955 series)

D

DIAGNOSTIC INSTRUMENT




L-2510HB

Synoptoscope with
Haidinger brush device



L-1550

Projection perimeter



KDM-3

"M-CHARTS"
Metamorphopsia
quantification charts



ISH-001

ISHIHARA colour
deficiency charts 38 plates

ISH-002: 24 plates
ISH-003: 10 plates
ISH-005: 14 plates



K-0161

Exophthalmometer
HERTEL

S

SURGICAL EQUIPMENT



S-5

Medical Quick Absorber



S-88ST

Medical Quick Absorber Stick Stout



G-701

Eye Shield Clear

DIAGNOSTIC INSTRUMENT



K-0391

Universal trial frame
Colour: K-0391 Black
K-0391-V Violet
K-0391-G White



K-350 (K-350A)

Trial lens set

35 PAIRS(Total 204pcs/228pcs) in wooden carrying case

SPHERES(-)	CYLINDERS(-)		CYLINDERS(+)		SPHERES(+)
35pairs	K-350 11pairs	K-350A 17pairs	K-350 11pairs	K-350A 17pairs	35pairs
20.00		6.00		6.00	20.00
16.00		5.00		5.00	16.00
13.00		4.50		4.50	13.00
12.00		4.00		4.00	12.00
11.00		3.50		3.50	11.00
10.00	3.00	3.00	3.00	3.00	10.00
9.00	2.50	2.50	2.50	2.50	9.00
8.00	2.25	2.25	2.25	2.25	8.00
7.00	2.00	2.00	2.00	2.00	7.00
6.50	1.75	1.75	1.75	1.75	6.50
6.00	1.50	1.50	1.50	1.50	6.00
5.75	1.25	1.25	1.25	1.25	5.75
5.50	1.00	1.00	1.00	1.00	5.50
5.25	0.75	0.75	0.75	0.75	5.25
5.00	0.50	0.50	0.50	0.50	5.00
4.75	0.25	0.25	0.25	0.25	4.75
4.50		0.12		0.12	4.50
4.25	AUXILIARIES (12 PIECES)		PRISMS (8 PIECES)		4.25
4.00					4.00
3.75	BLUE GLASS RED GLASS GREEN GLASS BROWN GLASS PLANE GLASS-2 OCCLUDER PINHOLE DISCS 0.5mm-1 1.0mm-1 1.5mm-1 SLIT DISCS 0.5mm-1 1.0mm-1		6 4 3 2 1 1 1/2 1/2		3.75
3.50					3.50
3.25					3.25
3.00					3.00
2.75					2.75
2.50					2.50
2.25					2.25
2.00					2.00
1.75					1.75
1.50					1.50
1.25					1.25
1.00					1.00
0.75					0.75
0.50					0.50
0.25					0.25
0.12					0.12

STAND



K-1411

Manual stand



K-1425-PT

Electric stand

Inami & Co., Ltd.

No.24-2, Hongo 3-chome, Bunkyo-ku,
Tokyo 113-0033, Japan
Tel: +81 (3) 3814-5916 Fax: +81 (3) 5684-2126
E-mail: corpinfo@inami.co.jp