

TOKYO

Inami

Cataract Surgery Set

www.inami.co.jp/en



Cataract Surgery Set



M-5R



M-1RW



S-511C



KI-46



M-901



M-1508



S-200S



M-1524



DS-50102



DS-1612F



M-167RA



DS-088



M-2177-5.5



M-136



M-136N



M-1160



M-1160N



DS-502



M-127C




M-128C

M-5R

Forceps Cornea and Suturing Type C Ultra Fine 0.12mm
length : 85mm, teeth width : 0.12mm




3D view



Ultra Fine

M-1RW

Forceps Suturing and Tying Wide handle
length : 87mm, tip width : 0.30mm length : 103mm





Wide Handle

S-511C

Scissors Spring Handle Type Curved
length : 105mm, blade length : 8mm



3D view



Utility

KI-46

Needle Holder for Microsurgery KOZAKI Curved without Lock
length : 105mm, tip width : 6mm




3D view



Dr. KOZAKI Model

M-901

Speculum for Microsurgery Open Ring model
lid length : 14mm



Standard

M-1508


Cannula Sub-Tenon OHMI Curved Three-Port 24G
length : 37mm, pipe : 24G, shaft length : 21mm



Tenon 3port

S-200S

Caliper CASTROVIEJO Straight Tip
length : 85mm



Standard

M-1524

Cannula Hydrodissection 23G
length : 40mm, pipe : 23G, shaft length : 29mm, teeth width : 0.7mm



45° Flat × 23G

DS-50102

Forceps Micro Capsulorhexis 23G Short Type
length : 136mm, shaft : 23G, shaft length : 15mm



CCC Side Port

DS-1612F


Forceps Capsulorhexis INAMURA Fine Titanium
length : 103mm



Dr. INAMURA Model

M-167RA

Forceps Circular Capsulorhexis UTRATA
length : 104mm



Dr. UTRATA Model

DS-088


Scissors YAE Multi-Purpose 23G 15mm
length : 137.5mm, shaft : 23G, blade length : 2.88mm



3D Action

M-2177-5.5

Marker CCC SUZUKI 5.5mm
length : 145mm, tip size : 5.5mm



CCC Marker

M-136

Hook Nucleus Dividing and Separating TOKUDA
length : 102mm, tip size : 0.6mm x 1.0mm



BUNKATSU-KUN

M-136N

PHACO CHOPPER NAGAHARA
length : 100mm, tip size : 0.3mm x 1.0mm



NAGAHARA

M-1160


Spatula Hook Nucleus Dividing and Separating Cobra-Shaft FUKUYAMA-YOSHITOMI
length : 104mm, tip size : 0.6mm x 1.0mm



COBRA-Spatula

M-1160N


Cobra-Shaft PHACO CHOPPER
length : 100mm, tip size : 0.3mm x 1.0mm



COBRA Chopper

DS-502

Micro Prechopper
length : 120mm, tip size : 1.5mm x 3.5mm



Prechopper

M-127C

Hook IOL SINSKEY Curved
length : 99mm, tip size : 0.2mm x 0.5mm



SINSKEY

M-128C

Hook IOL Curved Push and Pull Model
length : 109mm, tip size : 0.5mm x 0.7mm



Push & Pull



Craftsmanship

One of the cores of the skill that supports Inami's high quality products is "memory" of the craftsman's fingertips. Our craftsmen are not just routine workers but are always going through trials and errors for every improvement. And the results of the efforts are marked in their fingertips. The skills polished this way have been inherited generation to generation and support today's high quality of Inami products.

In order to maintain the "memory" and incarnate further quality improvement, we have also been concentrating our effort on training of young craftsmen and standardisation of the skills.



No.24-2, Hongo 3-chome, Bunkyo-ku,
Tokyo 113-0033, Japan
Tel: +81 (3) 3814-5916 Fax: +81 (3) 5684-2126
E-mail: corpinfo@inami.co.jp