

1. Class IIa or Im devices

EC REP
WMDE B.V.
Bergerweg 18
6085 AT Horn
The Netherlands



1-1. Tanito ab interno Trabeculotomy Hook Disposable (class IIa)

Model Name	Whole length	Shape	Body Color	Remarks
DP-2015S	126mm	Straight	Blue	Contains 5 pcs/box
DP-2015SET	126mm	Right/Straight/Left	Pink(Right)/Blue(Straight)/Yellow(Left)	5 sets of 3 pcs/box



1-2. Synoptoscope (generic name as Synoptophore) (class Im)

Model Name	Measurement	Weight	PD adjustment	Bulb	Heidenger Brush
L-2510B	550(W) x 320(D) x 392(H)	15kg	45 - 80mm	8V 0.75A	None
L-2510HB(with Heidenger Brush)	550(W) x 320(D) x 425(H)	15kg			50-100 round/min



1-3. Applanation Tonometer (class Im)

Model Name	Measurement	Diameter of the pneumatic face	Size	Weight	Standard Accessories
L-5110 Swing Arm Type	0 ~ 80mmHg in 2mm HG increments ±0.5mmHg	3.06mm	47mm (w)x 30mm (d)x 85mm (h)	725gs	Prism cone x 2, Caribration bar x 1, Adaptor x 1
L-5130 Plate Type					

2. Class I device

EC REP
WMDE B.V.
Bergerweg 18
6085 AT Horn
The Netherlands



Universal Light Weight Trial Frame

Model Name	Body Color	PD Adjustment	Cornea masurement	Frame angle	Frame Length	Max.number of lense
K-0391/K-0391T	Black	48-80mm	2mm interval	35°	35mm	10 pcs
K-0391G	White					
K-0391V	Violet					



Exophthalmometer, "HERTEL"

Model Name	Distance of precise protrusion measurement range	Distance of orbit-edge measurement range	Dimension	Weight
K-0161	0-35mm	78-112mm	76(D)x231(W)x19(H)mm	725g



Inami & Co., Ltd.

No.24-2, Hongo 3-chome, Bunkyo-ku,
Tokyo 113-0033, Japan
Tel: +81 (3) 3814-5916 Fax: +81 (3) 5684-2126
E-mail: corpinfo@inami.co.jp

Specifications are subject to change without prior notice.

Inami

Inami Ophthalmic Devices Instruments for Europe (EU)

