

K-0391/K-0391T

UNIVERSAL TRIAL FRAME

UNIVERSAL TRIAL FRAME K-0391/K-0391T



UNIVERSAL TRIAL FRAMES are designed to maximize examinee's comfort by minimizing frame weight yet retaining high durability. Quick and easy adjustment of both PD and corneal alignment are obtainable, while providing a sure mounting for trial lenses.

K-0391/K-0391T



Specifications:

PD : 48~80mm(± 0.8 mm)

Cornea measuring plate : 2mm interval

Frame adjustable angles : 35°(Upward 30° $\pm 5^\circ$, downward 5°+5°-2°)

Frame adjustable length : 35mm

Nose pad adjustable angles : 360°

Nose pad adjustable length : 22mm ± 1.5 mm

Maximum available lenses : 10pcs

Axis of lens : 0° $\pm 5^\circ$, 180° $\pm 5^\circ$

Available color : Black(K-0391/K-0391T), Violet(K-0391-V), White(K-0391-G)

A Passion for a Nice View



Inami & Co., Ltd.

No.24-2, Hongo 3-chome, Bunkyo-ku,
Tokyo 113-0033, Japan
Tel: +81 (3) 3814-5916 Fax: +81 (3) 5684-2126
E-mail: corpinfo@inami.co.jp

Specifications are subject to change without prior notice.



Management
System
EN ISO
13485:2016
www.tuv.com
ID 9000006746

EU REP
WMDE B.V.
Bergenvweg 18
6085 AT Horn
The Netherlands

