

NOGUCHI Lens Grabber II

Innovative ICL Extraction Forceps



Wins in Cataract complications/
Challenging cases !!




One-action IOL extraction Ideal for ICL setting

DS-2022S



DS-2022L



NOGUCHI Lens Grabber II

DS-2022S whole length 137mm /shaft 13.5mm /diameter 19G

DS-2022L whole length 159mm /shaft 36.0mm /diameter 19G

Material:Titun(tips:SUS)

Comment from Dr. Santaro Noguchi

(ASUCA Eye Clinic, Sendai City, Miyagi Prefecture/TSUKAZAKI Hospital, Ophthalmology,Himeji City, Hyogo Prefecture)

Intraocular lenses were basically extracted through a small incision after intraocular division.

However, with Lens Grabber, the entire lens can be removed in one action through a incision.

In addition, Lens Grabber was designed with a higher strength to increase the gripping force which was resulted in improving the one-action lens removal rate and decreasing forceps damage rates. Moreover, silicone lenses, which have been difficult to remove in the past, can now be removed in a single action without dividing the lens and in the long version, ICL can be easily and safely loaded into the cartridge like a Pac-Man forceps.




Hydrophilic Lens Extraction



Trifocal Lens Extraction



Conventional lens forceps



Lens Grabber II

