








Synoptiscope

L-2510B/HB

MADE IN JAPAN



Synoptiscope

L-2510B
L-2510HB

Inami Synoptiscope (generic name: Synoptophore) is designed for evaluations or trainings of binocular visions of strabismus or amblyopia with following three examinations:




S.P



Fusion



Stereopsis



L-2510B

Standard Version

L-2510HB

Haidinger Brush Version

Japan's precise machining and sturdiness of moulded body assure the long-lasting accuracy of examinations.



Detachable Haidinger Brush rotates at the speed of 50 to 100ppm.
It's motor has noise reduction structure.



Automatic flashing dials with 30 to 150 times per minutes are adopted along with manual flashing buttons on the control panel.



Side-painted pointer is easy and clear to read the scale.

For details see the website:



Inami & Co., Ltd.

No.24-2, Hongo 3-chome, Bunkyo-ku,
Tokyo 113-0033, Japan
Tel: +81 (3) 3814-5916 Fax: +81 (3) 5684-2126
E-mail: corpinfo@inami.co.jp



EU REP
WMDE B.V.
Bergerweg 18
6085 AT Horn
The Netherlands

CE
0197